Service- und Reparaturset

egger eLED.LP3

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REF 31700

REF 31710

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Wendemotor (REF 89320)
 Serien-Nr.:

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oder Wendemotor (REF 89330)
 Serien-Nr.:

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Service and Repair Kit egger eLED.LP3

5.7 Processing Times of the Individual LP/H Materials

LP/H Material	1st Polymerization	2nd Polymerization	3rd Polymerization
LP/H transparent	2 min. + black foil	3 min.	5 min. LP Cleaner
LP/H pink-transparent	2 min. + black foil	3 min.	5 min. LP Cleaner
LP/H beige	1 min. + black foil	3 min.	5 min. LP Cleaner
LP/H flesh-colored	1 min. + black foil	3 min.	5 min. LP Cleaner
LP/H blue	15 sec. + black foil	3 min.	5 min. LP Cleaner
LP/H red	15 sec. + black foil	3 min.	5 min. LP Cleaner

6 Warranty / Service / After-Sales Service

The warranty period for our devices starts with the date of delivery and is 12 months. Excluded from the warranty are modifications of the unit such as lamp modules, reflectors as well as all other parts of the unit which are subject to wear and tear.

In case of improper handling of the electric light-curing unit all warranty and liability claims will become extinct!

We will be pleased to be of assistance to you in case of damage or other malfunctions:

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E-mail international: sales@eggeronline.com

Internet: www.eggeronline.com

Please indicate the following data:

Device type:

Electric light-curing unit egger eLED.LP3 (REF 51350)

serial no.:

.....
 (see label)

Turning motor (REF 89320)

serial no.:

.....
 (see label)

or Turning motor (REF 89330)

serial no.:

.....
 (see label)

- Then, depending on the coating thickness applied, let it cure for about 2 – 6 minutes in the light-curing unit egger eLED.LP3.
- After the polymerization of the surface, wipe off the earmold by using a surface cleaner moistened cloth.
- If required, rework the earmold with a cutter and prepare it for lacquering as described in 5.3 / 5.4

5.6 Optional Manufacturing of BTE Earmolds and ITE Shells

Using the optionally available metalized insert, processing of egger LP/H-UV-polymerizates for the manufacturing of BTE earmolds and ITE shells is possible.

Tip: the metalized insert carton is fixed with included adhesive dots at insert. For this purpose, adhesive dots are attached to punched marking of metalized carton and pressed to the insert.

- Then, clean and degrease the surface carefully with surface cleaner by using the hand-spray pump (REF 86775).
 - Apply the lacquer evenly with the brush. Alternatively, let the lacquer dry in the light-curing unit rotating for about 6 – 9 minutes.
 - After the polymerization, the earmold (but do not immerse ITE earmolds with built-in electronics) is put into the LP-cleaner bath (in the round jar REF 52800) and then a second polymerization in the light-curing unit egger eLED.LP3 is to be done for about 6 minutes. This process is necessary to remove the inhibition layer (sticky surface) as a result from the first polymerization.
 - To finish the process, the earmold is rinsed under water and is blown dry.
 - Then, clean earmolds for 3 minutes in ultrasonic bath of water /isopropanol 1:1 (weighting %), treat with isopropanol and let dry at 37 °C for approx. 60 minutes.
- 5.5 Modifications / Repairs with Modeling Pastes**
- The LP/H modeling pastes can be universally used for various modifications and repairs with all kinds of hard UV-curing polymerizates or acrylates.

- First of all, the part of the earmold that is to be processed always has to be roughened with corundum paper, grit 180, (REF 40704).
- Then, clean and degrease the surface carefully with surface cleaner by using the hand-spray pump (REF 86775).
- Apply the modeling paste directly from the tube with the modeling instrument (REF 31215) and smooth the transitions to the ground.

5 Working with the UV-Polymerizates

5.1 Specific Storage Regulations

The materials are to be stored in a cool and dry place (see processing instructions for the various materials). Please mind that the containers are closed at any time, when using the LP/H lacquer antibac, you can protect hard earmolds with a special antibacterial coating with high gloss effect.

Otherwise, there might be danger of activating an unintentional chemical reaction through the permeability of light (curing of the material). The expiry date, applicable when stored adequately, is indicated on the product.

TIP: We recommend keeping already opened packages of lacquers or instant adhesives in the refrigerator.

5.2 Specific Processing Indications

Please avoid direct skin contact with the materials. Before using them protect your hands by using special skin protection products or wearing gloves.

The materials are to be used by authorized and trained staff only. Processing of the lacquers should be done underneath a suction unit. Always use the lacquering brushes for one type of lacquer only.

Details on the disposal for each product can be taken from our material safety data sheets, which can be downloaded in PDF format from our website:

www.egger-labor.com.

5.3 Lacquering with LP/H Lacquer – LP/H Lacquer antibac

For the surface finish all at present known hard acrylate-based otoplastic materials can be coated with LP/H lacquer.

Using the LP/H lacquer antibac, you can protect hard earmolds with a special antibacterial coating with high gloss effect.

The entire surface of the earmold always has to be roughened with corundum paper, grit 180, (REF 40704) first.

- Then, fix the BTE earmold onto a piece of nylon tube (REF 97200). For an ITE earmold, the alligator clip of the turning motor can be fixed on the opened battery door.
- Then, clean and degrease the surface carefully with surface cleaner by using the hand-spray pump (REF 86775).
- Apply the lacquer evenly with the brush. Alternatively, the earmold can also be immersed in a dipping jar (REF 31400 or 31401) filled with lacquer, but do not immerse ITE earmolds with built-in electronics or face-plates! Then, make sure that excess lacquer is shaken off.
- Let the lacquer dry in the light-curing unit rotating for about 2 – 3 minutes. Afterwards the earmold is directly ready for use.

5.4 Lacquering with LP/W Lacquer

For the surface finish all at present known hard acrylate-based otoplastic materials can be coated with LP/W lacquer to improve the fitting of the earmold ("non-skid coating").

- The entire surface of the earmold always has to be roughened with corundum paper, grit 180, (REF 40704) first.

- Position the device horizontally and place it in its suitable working position.
- Protect the unit against vibration.
- Check whether the unit is dry.
- Check whether the device is in technically perfect condition, especially the mains cable!
- Plug in the mains cable.

4.2 Initial Operating and Handling of the Light-Curing Unit egger eLED.LP3

Set the timer to the desired time (1, 3 or 5 minutes) by pressing the timer button and then start the polymerization process by pressing the start button.

Meaning of illuminating diodes

- 1 minute • • •
- 2 minutes • • • •
- 5 minutes • • • • •

Stopping the polymerization process

Keeping the start-stop-button pressed longer, the polymerization process can be stopped. In case of not using the device for 30 minutes, it changes into state of rest. The illuminating diodes turn off.

4.3 Advice for Disposal

The symbol on the unit means, that this unit may not be disposed in general trash, but must be disposed according to regional recycling guidelines for electronic devices and applicable environmental and disposal rules. For further information on disposal and recycling please contact your regional authorities or your distributor of the unit.

- Use the unit in perfect and clean condition only and by observing the operating instructions.
- Connect the device to permissible sockets only.
- Use the included cable only.
- Please adhere to the acknowledged safety regulations.
- This device is used for the processing of UV-curing materials in a spectral range of 360 – 400 nm.
- Please also adhere to the processing instructions for each material used.
- Disconnect the unit from electricity supply prior to any service or maintenance.
- Do not pull the mains plug out of the socket at the mains cable.
- Should the mains cable be damaged, it must be replaced at once. Defective mains cables must not be used.
- The unit may not be opened by operators. There are no interchangeable parts such as fuses inside.
- Repairs may only be executed by the manufacturer.

Advice for cleaning: Clean the unit with a moist cloth only – do not use solvents, do not douse the unit into water.

On self-accomplished repairs or if other spare parts are used, all warranty and liability claims will become extinct!

4 Installation, Handling, Initial Operation

4.1 Assembling and Connecting the Unit

The electrical light-curing unit egger eLED.LP3 corresponds to the state-of-the-art standards and to the acknowledged safety regulations. Nevertheless dangers can arise. Please observe the following, before putting the unit into operation:

- The device is to be positioned in a way that a fast separation of mains is possible very easily and at any time.

- REF 29615 1 tube LP/H modeling paste blue 20 g
- REF 306230 incl. 5 disposable brushes
- REF 306230 1 bottle LP/H lacquer antibac 30 ml,
- REF 307030 incl. 5 disposable brushes
- REF 31301 1 bottle LP-cleaner 100 ml
- REF 31600 1 bottle surface cleaner 250 ml
- REF 86775 1 hand-spray pump for surface cleaner
- REF 52800 1 jar round, for LP-cleaner
- REF 97200 1 meter nylon tube as support for lacquering
- REF 31215 1 modelling instrument
- REF 68601 1 timer to control polymerization times
- REF 39503 1 Uniglu®/2 instant adhesive
- REF 39504 5 measuring tips large, for Uniglu®
- REF 39505 5 measuring tips small, for Uniglu®
- REF 95100 10 sound tubes pre-bent, 2 x 3.1 mm
- REF 95100 1 user and processing instruction

The electrical light-curing unit egger eLED.LP3 is used for the polymerization of UV-curing modelling pastes and lacquers. A use other than described in the operating instructions may lead to injuries of the operator or of third persons.

Furthermore, the device or other material assets may be damaged. The device must be used in technically perfect condition only and the safety regulations must be strictly complied with.

Warning hot surfaces – do not touch the LED reflectors – risk of burns!

Work area identification, Risk class 3, according to DIN EN 62471.

At a distance of less than 20 cm from the lamp, the exposure limits for UV radiation are exceeded.

For the egger eLED.LP3 light polymerization unit, the operator must prepare a hazardous situation assessment and safety instructions for operation.

2.2 Technical Data

Power supply: U = 100 – 240 V (AC)
 Mains frequency: f = 50/60 Hz
 Power input: P = 22 VA
 Life expectancy of the LED bulbs: 20,000 operating hours
 Connection: Euro plug

3.1 For a Safety Conscious Work Environment

- Protect the unit from dust. Use the unit on dry working places only.
- Store the unit in closed rooms.
- Qualified and trained persons must use the device only.
- Children must be kept away from the unit.
- Surrounding temperature between 10 °C and 30 °C.
- Relative air humidity between 10 % and 70 %.
- Before putting the device into operation, the user must read the operating instructions and ask the manufacturer in case of any questions.

The light-curing unit egger eLED.LP3 has been especially developed for the above mentioned applications and is characterized by

- Maintenance-free working
- Durable LED light sources
- Easy handling
- Simple cleaning and care of the powder-coated metal casing

This egger service and repair kit eLED.LP3 is rounded off by the necessary accessories for a quick gluing in of sound tubes.

Furthermore, we can offer a variety of devices and useful accessories from our wide range of products such as:

- modeling pastes, lacquers and accessories
- cutters, drills und tools
- technical units
- compressors for compressed air
- compressed-air sprays

2 Scope of delivery, Accessories, Technical Data

2.1 Scope of Delivery

Please check whether the delivery is complete and whether all parts are free from transport damage. Immediately inform us or the distributor about missing parts or transport damage.

The service and repair kit eLED.LP3 consists of:

- REF 51350 1 light unit egger eLED.LP3
- REF 89320 1 turning motor battery-operated (kit 31700) REF 89320 or 1 electrical turning motor (kit 31710) REF 89330
- REF 30800 1 tube LP/H modeling paste transparent 20 g REF 30800
- REF 31000 1 tube LP/H modeling paste beige 20 g REF 31000
- REF 29610 1 tube LP/H modeling paste red 20 g REF 29610

1 General

The egger service- and repair kit is suited for a number of uses in the acoustician shop, during home visits and in the optoelastic laboratory. Service operations can be executed immediately right in the shop – fast, efficiently, no lengthy preparation needed. So you can enjoy more independence while offering your customer the best possible service.

All possible uses and advantages at one glance:

- ### Modeling pastes
- Repairs of BTE and ITE earmolds that are made of UV-polymerizates or other hard acrylicates
 - Minimizing or sealing ventilation bores
 - Post-production attachment of retention points or a lip seal
 - Modifications of vents on ITE earmolds
 - Handy 20 g tubes for efficient and clean processing
- ### Lacquers
- Coatings with LP/H lacquer or LP/H lacquer antibac
 - Post-production varnishing ("glazing") of polished earmolds using LP/H lacquers – a preventive measure against skin irritation
 - Partial build-up of earmolds to improve their seal
 - Adhesive coating ("non-slip coating") using LP/W lacquer to improve an earmold's hold
 - Usable for all UV-polymerizates and hard acrylicates
 - Available in handy 30 ml bottles including lacquering brush

The service and repair kit eLED.LP3 consists of:

- REF 51350 1 light unit egger eLED.LP3
- REF 89320 1 turning motor battery-operated (kit 31700) REF 89320 or 1 electrical turning motor (kit 31710) REF 89330
- REF 30800 1 tube LP/H modeling paste transparent 20 g REF 30800
- REF 31000 1 tube LP/H modeling paste beige 20 g REF 31000
- REF 29610 1 tube LP/H modeling paste red 20 g REF 29610

Please read the operating instructions carefully before working with the device!
Please keep the operating instructions in a safe place for future reference!

Index

1 General

2 Scope of delivery, Accessories, Technical Data

3 Safety Regulations

4 Installation, Handling, Initial Operation

5 Working with UV-Polymerizates

6 Warranty / Service / After-Sales Service

7

8

Caution

Warning hot surfaces – risk of burns!

Work area identification, Risk class 3, according to DIN EN 62471. Caution dangerous ultraviolet radiation!

Class II equipment.

Caution!

Important Information!

Warning hot surfaces – risk of burns!

Work area identification, Risk class 3, according to DIN EN 62471. Caution dangerous ultraviolet radiation!

Class II equipment.

Caution!

9

10

11

12

13

14

15

16

17

18

19

20

21

22

23

24

25

26

27

28

29

30

31

32

33

34

35

36

37

38

39

40

41

42

43

44

45

46

47

48

49

50

51

52

53

54

55

56

57

58

59

60

Service and Repair Kit

egger eLED.LP3

Operating instructions

REF 31700 (100 – 240 V)
REF 31710 (100 – 240 V)

egger

egger Otoplastik + Labortechnik GmbH
Aybühlweg 59
87439 Kempen / Germany