egger Acryl/B

Kaltpolymerisat
Cold curing acrylic resin
Šaltieji polimerai

egger

1. Hersteller

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Dentaurum GmbH & Co. KG
Turnstr. 31 · 75228 Ispringen · Germany

2. General device description

Acryl/B are polymethyl methacrylate (PMMA)-based cold-curing resins in various colors for the production of BTE earmolds and ITE shells. Direct use on patients is not permitted.

3. Intended purpose

Acrylics for the production of earmolds.

4. Intended users

Acryl/B may only be used by specialized earmold laboratories or by trained hearing care professionals.

5. Indications

Earmolds are components for hearing systems.

6. Contraindications

The fabricated earmold products must not be used for patients with known allergies to one of the product constituents (PMMA).

7. Production of a BTE earmold or ITE shell

7.1. Mixing

Mix the monomer (fluid) and polymer (powder) in a weight ratio of 1:2. The monomer should be cold (refrigerator temperature, 2 to 8 °C) and should be added to the mixing cup first when mixing. The mixture should be prepared on a vibrating plate to avoid inclusion of bubbles. After mixing, the acrylic must be processed immediately.

7.2. Fabrication of BTE earmolds

Pour the mixed and bubble-free acrylate mixture into the negative mold of the ear impression.

7.3 Fabrication of ITE hollow shells

The negative mold is heated to 60 °C. Pour the mixed and bubble-free acrylate mixture into the preheated negative mold. A shell wall forms after a few seconds. Then turn the mold upside down to allow the excess acrylic mixture to run off.

7.4. Polymerization

As soon as a matt surface has formed, the earmold can be polymerized. Heat the mold under water (40 °C) at 2–5 bar for 12–15 minutes.

7.5. Post treatment

If required, the shape and surface can then be reworked using suitable milling and polishing agents.

7.6. Varnishing

The earmold can then be coated with a CE-marked lacquer suitable for this intended purpose (e.g. egger LP/H lacquer).

8. Safety information

H225-H315-H317-H335

Highly flammable liquid and vapor. Causes irritation to the skin. Can cause allergic reactions. Can cause irritation to the respiratory tracts.

P210-P261-P280-P302+P352-P333+P313-P403

Keep away from heat, hot surfaces, sparks, open flames as well as other ignition sources. Do not smoke. Avoid breathing dust / fume / gas / mist / vapors / spray. Wear protective gloves / protective clothing / eye protection / face protection / hearing protection. IF ON SKIN: Wash with plenty of water. If skin irritation or rash occurs: Get medical advice/attention. Store in a well-ventilated place.

8. Information for single use products

The product is intended for single use only. Reconditioning (recycling) of the product and its reuse is not permitted, as perfect functionality of the product cannot be guaranteed.

9. Additional information

Should the user and/or the patient become aware of serious problems arising from the use of the product, it is important that the manufacturer and the competent authority in the country in which the user and/or the patient is resident are informed accordingly.

10. Notes on disposal

When disposing of the product, please adhere to the national regulations that are valid for your region and – if available – the details outlined in the safety data sheets.

11. Explanation of symbols used on the label

Name and address of the manufacturer:
 DENTAURUM GmbH & Co. KG
 Turnstr. 31 · 75228 Ispringen / Germany

Name and address of the distributor:
 egger Otoplastik + Labortechnik GmbH
 Aybühlweg 59 · 87439 Kempten / Germany

Protect against sunlight

Expiry date

HIBC symbol, for the structure
 of the barcode according to HIBC

0483 Manufacturer's CE mark of conformity

Refer to Instructions for use

Medical device

UDI symbol, for the carrier
 of the Unique Device Identifier

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