Sicherheitshinweise / Safety regulations:

DE: Achtung

EN: Warning

H: 315 - 319 - 335

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EN: Causes skin irritation. Causes serious eye irritation. May cause respiratory irritation.

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EN: Keep out of reach of children. Avoid breathing mist/vapours/spray. Use only outdoors or in a well-ventilated area. Wear protective gloves/protective clothing/eye protection. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Call a POISON CENTER/doctor if you feel unwell. Store locked up. Dispose of contents/container to hazardous or special waste collection point.

EUH202 Cyanoacrylate. Danger. Bonds skin and eyes in seconds. Keep out of the reach of children.

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Please follow the product information and the material safety data sheet.

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General product information

unigluе/2 is a rapid-hardening, solvent-free single component adhesive. This instant adhesive is used to bond different kinds of earmold materials. Please note that the surfaces of silicone earmolds should be pre-treated with primer/2 (REF 39701) before gluing sound tubes.

EUH202! Cyanoacrylate. Danger. Bonds skin and eyelids in seconds. Keep out of reach of children!
Contains: Ethyl 2-cyanoacrylate.

Storage and durability:

The shelf life is, at optimum storage temperature of +2 °C to +10 °C in the closed original container, a maximum of 24 months. A higher storage temperature leads to a significantly lower durability. The storage temperature must not fall below +2 °C. Protect the product from direct sunlight or heat, frost and humidity.

Usage:

1. Close container properly after usage.
2. Avoid contact with skin. Before handling, protect hands with appropriate gloves.

Handling:

Apply unigluе/2 always very sparingly. The thinner the bonded joint, the firmer the bond and the more quickly hardening takes place. The ultimate strength of the bond is achieved within 12 hours.

The usage of unigluе/2 requires handling by qualified personnel considering the relevant safety data sheet. Working with unigluе/2 should be done under a vent. Please see safety data sheet of this product for disposal.

Handling instruction:

1. Preparation

Best results are achieved on clean surfaces. Remove dust and harder particles by using compressed air. It is advisable to clean the surface beforehand with Surface Cleaner (REF 31600 ff) from greasy and oily remainders. Please consider the evaporation times of the used cleaner. For the generation of compressed air we offer various compressors and our cedis AirPower spray. The adhesive and the substrates to be bonded should be at a temperature between +18 °C and +25 °C before processing to optimize the adhesive strength.

2. Gluing

Apply a thin film of unigluе/2. Depending on the material and the amount of glue used, press the parts to be glued together tightly. Due to the capillary action the adhesive penetrates the smallest gaps. Advice: For the application of the unigluе/2 instant adhesive we offer dosing tips in the sizes 0.5 mm (REF 39505) and 1.0 mm (REF 39504).