**e3D.M01**

Otoplastiken & IdO Schalen  
Earmoulds & ITE Shells

**egger**

<b>DE</b>	Gebrauchsanweisung .....	6	<b>PT</b>	Manual de instruções .....	45
<b>EN</b>	Instructions for use.....	9	<b>RO</b>	Instrucțiuni de utilizare.....	48
<b>FR</b>	Mode d'emploi.....	12	<b>SV</b>	Bruksanvisning .....	51
<b>ES</b>	Modo de empleo .....	15	<b>SK</b>	Návod na použitie.....	54
<b>IT</b>	Istruzioni per l'uso.....	18	<b>SL</b>	Navodila za uporabo.....	57
<b>DA</b>	Brugsanvisning .....	21	<b>CS</b>	Návod k použití .....	60
<b>EL</b>	Οδηγίες χρήσης.....	24	<b>HU</b>	Használati útmutató.....	63
<b>HR</b>	Upute za uporabu.....	27	<b>RU</b>	Инструкция по применению .....	66
<b>LV</b>	Lietošanas instrukcija .....	30	<b>TR</b>	Kullanım kılavuzu .....	69
<b>LT</b>	Naudojimo instrukcija .....	33	<b>FI</b>	Käyttöohje .....	72
<b>NL</b>	Gebruiksaanwijzing .....	36	<b>ET</b>	Kasutusjuhend .....	75
<b>NO</b>	Bruksanvisning .....	39	<b>BG</b>	Инструкция за употреба .....	78
<b>PL</b>	Instrukcja stosowania .....	42			

## Suitable for the following DLP printers / cleaning / post-exposure

Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
<b>Asiga Max</b>	385 nm	6.1 mW/cm <sup>2</sup>	Composer 1.2.11	e3D.M01_Color
<b>Asiga Pico 2</b>	385 nm	20 mW/cm <sup>2</sup>	Composer 1.2.11	e3D.M01_Color
<b>Asiga PRO 2</b>	385 nm	5.7 mW/cm <sup>2</sup>	Composer 1.2.11	e3D.M01_Color
<b>Asiga PRO 4K</b>	385 nm	7.0 mW/cm <sup>2</sup>	Composer 1.2.11	e3D.M01_Color
<b>Microlay Versus</b>	385 nm	4.3 mW/cm <sup>2</sup>	Microform 1.0.3.7	e3D.M01_Color
<b>Miicraft Ultra Series</b>	385 nm	5.7 mW/cm <sup>2</sup>	Utility 6.3.0	e3D.M01_Color
<b>Way2Production SolFlex Series</b>	385 nm	8.0 mW/cm <sup>2</sup>	Netfabb 2020	e3D.M01_Color

\* The set of parameters includes all relevant material- and printer specific information.

The placeholder COLOR shall be replaced by **transparent, pink-transparent, neon-yellow, red, blue, blue plus** or **light-brown**

## Cleaning equipment

Cleaning unit Manufacturer, Model	Cleaning process
Ultrasonic bath Bandelin Sonorex	<p>Clean the parts with isopropyl alcohol (purity <math>\geq</math> 98 %) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air.</p> <p>The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity <math>\geq</math> 98 %) for 3 minutes.</p> <p>Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.</p>

## Light curing equipment

Light Curing unit Manufacturer, Model	Curing process
NK Optik Otoflash G171	2 x 2000 flashes under inert gas, turnaround components after 2000 flashes.

## Mechanical finishing

Polishing unit Manufacturer, Model	Polishing process
Otec ECO-Maxi wet	<p><b>Step 1:</b> Media: DZS 10/10 Ceramic grinding wheel, Compound: SC 15 Compound, Runtime: 90 min, Splitting system: Wet splitting, Speed: 280 U/min, Compound content: 3 %, Water flow: 4 l/h</p> <p><b>Step 2:</b> Media: KM 10 plastic abrasive media (70 %), PM 10 plastic abrasive media (30 %), Compound: SC 15 Compound, Runtime: 90 min, Splitting system: Wet splitting, Speed: 280 U/min, Compound content: 3 %, Water flow: 4 l/h</p> <p><b>Step 3:</b> Media: ZSP 4/5 Porcelain polishing body, Runtime: 60 min, Splitting system: Wet splitting, Speed: 260 U/min, Water flow: 10 l/h</p>
Otec CF 3x18	<p><b>Step 1:</b> Media: DZS 10/10 Ceramic grinding wheel, Compound: SC 15 Compound, Runtime: 90 min, Rib shape: standard, Splitting system: Wet splitting, Speed: 280 U/min, Compound content: 3 %, Water flow: 8 l/h</p> <p><b>Step 2:</b> Media: KM 10 plastic abrasive media (70 %), PM 10 plastic abrasive media (30 %), Compound: SC 15 Compound, Runtime: 90 min, Rib shape: standard, Splitting system: Wet splitting, Speed: 280 U/min, Compound content: 3 %, Water flow: 8 l/h</p> <p><b>Step 3:</b> Media: ZSP 4/5 Porcelain polishing body, Runtime: 60 min, Rib shape: standard, Splitting system: Wet splitting, Speed: 240 U/min, Water flow: 15 l/h</p>

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## Intended use

Resin material for the production of earmoulds and ITE shells in 3D printing processes

## Indication

Earmoulds for ear defenders and hearing aid earmoulds for connection to a hearing system, In-Ear-Monitoring

## Patient Target Group

Persons for whom an earmould is to be created.

## Intended Users

Earmould laboratory, ENT specialist, hearing aid acoustician

## Suitable for the following DLP printers / cleaning / post-exposure / Polishing machines

See page 4/5

## Processing

- The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the earmoulds are completely cured (refer to the process description).
- Bottle: Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- Eco Bag: Before use, the material should be shaken intensively and homogenized by tumbling or using a roller mixer with a suitable attachment.
- The final colour shade is achieved after a few days.
- Processing temperature  $23\text{ °C} \pm 2\text{ °C}$  /  $73,4\text{ °F} \pm 35,6\text{ °F}$ .
- During the first days after their fabrication don't expose the finished earmoulds to the direct sunlight.

## Safety Information

- Only for the specified intended use by trained specialists.
- Avoid direct contact with the liquid material and the components before post-curing, in particular in pregnant / breastfeeding women. Irritating to eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- In case of contact with skin, wash immediately with plenty of soap and water.
- Biocompatibility is guaranteed only after complete polymerization.
- The processing instructions and precautionary measures must be observed.
- Refer to the relevant safety data sheet for hazard and safety information.

## Important notes

- DETAX shall not be held liable for any damage caused by misuse.
- Always keep container tightly sealed, immediately close the container carefully after each use.
- Read and understand the safety data sheet!

### For users and /or patients:

Any serious incidents occurring in relation to this product should be reported immediately to [incident@detax.com](mailto:incident@detax.com) and to the competent authority of the Member State in which the user and/or patient is established.

## Storage

- e3D.M01 is to be stored dry (at 15–28 °C / 59–82 °F) and protected from light. Even slight exposure to light can trigger polymerization.
- If the material is stored at too low temperatures, it may become cloudy, but this can be reversed without loss of quality by warming the material briefly to 35 °C / 95 °F.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

## Contraindication

Contains (meth)acrylics and phosphine oxides.

Some ingredients of e3D.M01 may cause allergic reactions in predisposed persons. In such cases refrain from using the product.

## Adverse effects

Product may cause allergic reactions.

## Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

### Storage:

### Application:

At 23 °C ± 2 °C / 73 °F ± 4 °F

## Manufacturing process for earmoulds

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

### Construction process

Generation of a Print Job complying with machine and material parameters

### Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended.

If possible, post-processing should commence immediately following the construction process.

### Cleaning

See page 5

### Drying

Heat the construction components for 30 min. in an oven to approx. 40 °C / 104 °F to remove the solvent residues from the cleaning process.

### Post-exposure

See page 5

### Surface processing

#### Option 1

Varnish the surface with

LP/H Lacquer

#### Option 2

Mechanical finishing of the surface

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egger Otoplastik+ Labortechnik GmbH
DETAX GmbH