eWA.X3

Service-Station

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eWA.X3 Service Station

Developed in cooperation with hearing aid acousticians, the eWA.X3 is indispensable after only a short time when deposits (i. e. cerumen) or moisture are difficult to access and when they block sound's way. The device is useful for many applications: from receiver and microphone openings to sound tubes and slim tubes to cerumen protection sieves to additional bore holes, etc. To ensure safe, effective and gentle application, the integrated suction and air-blowing function has undergone intensive testing and has been perfectly matched to the everyday life at the service workplace. The perceptibly low noise level is the result of a careful selection of components, thought-out construction and a particular installation of the high-performance pump.

Initial Operation

In order to start up the eWA.X3, plug it in an electrical outlet. The device can be turned on and off via the rocker switch. When the device is turned on, the manometer should display about 0.1 bar when the tube ends are open (e.g. the attachment has been removed).

The value changes depending on the attachment.

For suction, choose the suitable attachment and plug it on top of the dirt particle filter. Now the cleaning process can commence. Treat soiled parts similarly to using a vacuum cleaner.

IMPORTANT:

Perform work using the service station cautiously in order to avoid damaging the hearing system or the earmold with the attachments.

Cleaning

When the value on the display rises to 0.7 bar or more, the attachment and/or the particle filter must be cleaned; it could also be the case that the suction tube is blocked.

The following elements are to be cleaned regularly in an ultrasonic bath:

- Attachments – to be cleaned after every use
- Dirt particle filter – depending on visual state, but to be cleaned at least once per week


In order to do so, pull the particle filter off the suction tube. Unscrew the transparent protection cover of the filter and pull out the filter. Put all of it into an ultrasonic bath. Take care of the sealing rings. When reassembling, they have to be reinserted correctly. Push the compounded filter – black part first – into the tube.

- Tubes – depending on visual state, but to be cleaned at least once per month

In order to clean the tubes, pull them off the suction device and put them into an ultrasonic bath. When fitting them together again, take care that there is no water left inside the tubes.

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Elements of the Device

1. Manometer
2. On/Off switch
3. Outgoing air tube ⇨
4. Suction tube ⇨
5. Dirt particle filter at the suction tube
6. Power plug at the reverse side of the device

Attachments

- purple Ø 0.51 mm
- pink Ø 0.58 mm
- green Ø 0.84 mm
- orange Ø 1.37 mm

Technical Data

- Voltage input:
AC 100–240 V; 50–60 Hz; 0.7 A
- Power consumption: 12 W
- Fuse: 2 A (sluggish)
- Dimensions: approx. W 178 x D 153 x H 70 mm
- Weight: 2.5 kg
- Volume: 12–15 LPM (liter per minute)
- Vacuum: -70 kPa
- Overpressure: 150 kPa
- Temperature range: 0–50 °C / 0–122 °F

Scope of Delivery

- eWA.X3 (REF 62850)
- 1 set of tubes
- 1 set of attachments (purple, pink, green, orange)
- 1 sintered filter

Spare Parts

- Set of tubes (REF 62870)
- Set of attachments (REF 62871)
- Sintered filter (REF 62872)