e3D.S01

IdO Schale
ITE Shells

egger

DE	Gebrauchsanweisung	6	PT	Manual de instruções	45
EN	Instructions for use.....	9	RO	Instrucțiuni de utilizare.....	48
FR	Mode d'emploi.....	12	SV	Bruksanvisning	51
ES	Modo de empleo	15	SK	Návod na použitie.....	54
IT	Istruzioni per l'uso.....	18	SL	Navodila za uporabo.....	57
DA	Brugsanvisning	21	CS	Návod k použití	60
EL	Οδηγίες χρήσης.....	24	HU	Használati útmutató.....	63
HR	Upute za uporabu.....	27	RU	Инструкция по применению	66
LV	Lietošanas instrukcija	30	TR	Kullanım kılavuzu	69
LT	Naudojimo instrukcija	33	FI	Käyttöohje	72
NL	Gebruiksaanwijzing	36	ET	Kasutusjuhend	75
NO	Bruksanvisning	39	BG	Инструкция за употреба	78
PL	Instrukcja stosowania	42			

Suitable for the following DLP printers / cleaning / post-exposure

Printer Manufacturer, Model	Light source	Light intensity	Operation Software	Parameter data set*
Asiga Max	385 nm	6.1 mW/cm ²	Composer 1.2.11	e3D.S01_Color
Asiga Pico 2	385 nm	20 mW/cm ²	Composer 1.2.11	e3D.S01_Color
Asiga PRO 2	385 nm	5.7 mW/cm ²	Composer 1.2.11	e3D.S01_Color
Asiga PRO 4K	385 nm	7.0 mW/cm ²	Composer 1.2.11	e3D.S01_Color
Microlay Versus	385 nm	4.3 mW/cm ²	Microform 1.0.3.7	e3D.S01_Color
Miicraft Ultra Series	385 nm	5.7 mW/cm ²	Utility 6.3.0	e3D.S01_Color
Way2Production SolFlex Series	385 nm	8.0 mW/cm ²	Netfabb 2020	e3D.S01_Color

* The set of parameters includes all relevant material- and printer specific information.
The placeholder COLOR shall be replaced by **white**, **beige**, **black**, **red plus** or **blue plus**.

Cleaning equipment

Cleaning unit Manufacturer, Model	Cleaning process
Ultrasonic bath Bandelin Sonorex	<p>Clean the parts with isopropyl alcohol (purity $\geq 98\%$) for 3 minutes. Then thoroughly clean the openings, cavities and gap areas with compressed air.</p> <p>The main cleaning is performed in a separate vessel with fresh isopropyl alcohol (purity $\geq 98\%$) for 3 minutes.</p> <p>Prior to post-exposure, check the openings, cavities and gap areas for residues. Then blow off with compressed air.</p>

Light curing equipment

Light Curing unit Manufacturer, Model	Curing process
NK Optik Otofash G171	2 x 2000 flashes under inert gas, turnaround components after 2000 flashes (expose the color black with 3 x 2000 flashes)

Mechanical finishing

Polishing unit Manufacturer, Model	Polishing process
Otec ECO-Maxi wet	<p>Step 1: Media: DZS 10/10 Ceramic grinding wheel, Compound: SC 15 Compound, Runtime: 90 min, Splitting system: Wet splitting, Speed: 280 U/min, Compound content: 3 %, Water flow: 4 l/h</p> <p>Step 2: Media: KM 10 plastic abrasive media (70 %), PM 10 plastic abrasive media (30 %), Compound: SC 15 Compound, Runtime: 90 min, Splitting system: Wet splitting, Speed: 280 U/min, Compound content: 3 %, Water flow: 4 l/h</p> <p>Step 3: Media: ZSP 4/5 Porcelain polishing body, Runtime: 60 min, Splitting system: Wet splitting, Speed: 260 U/min, Water flow: 10 l/h</p>
Otec CF 3x18	<p>Step 1: Media: DZS 10/10 Ceramic grinding wheel, Compound: SC 15 Compound, Runtime: 90 min, Rib shape: standard, Splitting system: Wet splitting, Speed: 280 U/min, Compound content: 3 %, Water flow: 8 l/h</p> <p>Step 2: Media: KM 10 plastic abrasive media (70 %), PM 10 plastic abrasive media (30 %), Compound: SC 15 Compound, Runtime: 90 min, Rib shape: standard, Splitting system: Wet splitting, Speed: 280 U/min, Compound content: 3 %, Water flow: 8 l/h</p> <p>Step 3: Media: ZSP 4/5 Porcelain polishing body, Runtime: 60 min, Rib shape: standard, Splitting system: Wet splitting, Speed: 240 U/min, Water flow: 15 l/h</p>

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Intended use

Resin for the production of ITE shells in 3D printing processes

Indication

In-the-ear shells for connection to a hearing system

Patient Target Group

Persons for whom an earmould is to be created.

Intended Users

Earmould laboratory, ENT specialist, hearing aid acoustician

Suitable for the following DLP printers / cleaning / post-exposure / Polishing machines

See page 4/5

Processing

- The properties of the final product depend, among other things, on post-processing. Correct post-exposure is important for biocompatibility. Therefore it must be ensured that the light unit is in an orderly condition and that the earmoulds are completely cured (refer to the process description).
- Before use, the material in the bottle should be shaken intensively and homogenized with a bottle roller.
- The final colour shade is achieved after a few days.
- Processing temperature $23\text{ °C} \pm 2\text{ °C} / 73,4\text{ °F} \pm 35,6\text{ °F}$.

Safety Information

- Only for the specified intended use by trained specialists.
- Avoid direct contact with the liquid material and the components before post-curing, in particular in pregnant/breastfeeding women. Irritating to eyes and skin (sensitization possible).
- Wear personal protective equipment (protective gloves, goggles) when handling the uncured material.
- Wear suitable personal protective equipment (protective gloves, goggles, face mask) when finishing the cured material.
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.
- In case of contact with skin, wash immediately with plenty of soap and water.
- Biocompatibility is guaranteed only after complete polymerization.
- The processing instructions and precautionary measures must be observed.
- Refer to the relevant safety data sheet for hazard and safety information.

Important notes

- DETAX shall not be held liable for any damage caused by misuse.
- Always keep container tightly sealed, immediately close the container carefully after each use.
- Read and understand the safety data sheet!

For users and /or patients:

Any serious incidents occurring in relation to this product should be reported immediately to incident@detax.com and to the competent authority of the Member State in which the user and/or patient is established.

Storage

- e3D.S01 is to be stored dry (at 15–28 °C / 59–82 °F) and protected from light. Even slight exposure to light can trigger polymerization.
- To protect against contamination, cover the material in the tray with the lid or a glass plate.

Contraindication

Contains (meth)acrylics and phosphine oxides.

Some ingredients of e3D.S01 may cause allergic reactions in predisposed persons. In such cases refrain from using the product.

Adverse effects

Product may cause allergic reactions.

Disposal

Disposal of the contents /container must be carried out in accordance with the local /regional /national and international regulations.

Storage:

Application:

At 23 °C ± 2 °C / 73 °F ± 4 °F

Manufacturing process for earmoulds

Data preparation and fabrication of the support structure according to the instructions of the CAD software manufacturer

Construction process

Generation of a Print Job complying with machine and material parameters

Post-processing

After raising the platform, a drip time of approx. 10 minutes is recommended.

If possible, post-processing should commence immediately following the construction process.

Cleaning

See page 5

Drying

Heat the construction components for 30 min. in an oven to approx. 40 °C / 104 °F to remove the solvent residues from the cleaning process.

Post-exposure

See page 5

Surface processing

Option 1

Varnish the surface with
LP/H Lacquer

Option 2

Mechanical finishing of the surface
(polishing, vibratory finishing)

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egger Otoplastik+ Labortechnik GmbH
DETAX GmbH