**egger A vario**  
**egger A/I vario**  
**egger A/II vario**

Additionsvernetzendes  
Ohrabformsilikon

Addition curing ear  
impression silicone

**egger**

## Technical Specifications

|                                                   | egger A vario                                                                     | egger A/I vario                                                                   | egger A/II vario                                                                   |
|---------------------------------------------------|-----------------------------------------------------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|
| <b>Product color</b><br>(Base/catalyst)           | blue/white                                                                        | green/white                                                                       | red/white                                                                          |
| <b>Mixing volume</b>                              | 2 x 700 g                                                                         | 2 x 660 g                                                                         | 2 x 700 g                                                                          |
| <b>Texture, consistency</b><br>low high           |  |  |  |
| <b>Open time</b><br>Min. approx.                  | 2:00                                                                              | 2:00                                                                              | 2:00                                                                               |
| <b>Indwelling time in the ear</b><br>Min. approx. | 3:00                                                                              | 3:00                                                                              | 3:00                                                                               |
| <b>Total setting time</b><br>Min. approx.         | 5:00                                                                              | 5:00                                                                              | 5:00                                                                               |
| <b>Final hardness</b>                             | 35 Shore A                                                                        | 30 Shore A                                                                        | 40 Shore A                                                                         |

The processing time refers to a room temperature of  $23^{\circ}\text{C} \pm 1^{\circ}\text{C}$  ( $73.4^{\circ}\text{F} \pm 1.8^{\circ}\text{F}$ ) and a relative air humidity of 50 %, the setting time to a body temperature of  $37^{\circ}\text{C} \pm 1^{\circ}\text{C}$  ( $98.6^{\circ}\text{F} \pm 1.8^{\circ}\text{F}$ ). Technical data according to DIN EN ISO 4823.

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**Caution!**

The impression of the outer ear canal may only be taken with impression material by persons with appropriate qualifications. The processing information and precautions which are set out in these directions for use must be observed under all circumstances. Failure to observe these may result in irreparable damage to the ear or the eardrum. If damage or undesirable side effects occur despite the warnings, we request that you report this with a detailed description of your procedure and the damage to or side effects for the patient.

**Intended use**

Addition curing ear impression silicone

**Indications for use**

Ear impression material for reproduction of the outer ear canal of a patient with hearing loss, Tinnitus aurium or impaired hearing for the provision and adjustment of a hearing device.

**Contraindication**

The product may not be used if the patient is allergic to one of the ingredients. There are no side-effects to be expected from this medical device if processed and used properly. However, immune reactions (e.g. allergies) or local disorders of sensation (e.g. irritation in the auditory canal) cannot be ruled out in principle. If you become aware of adverse side-effects – and in case of doubt – please be sure to advise us in any case and describe the circumstances and symptoms as accurately as possible. We will investigate any indication.

**Patient groups**

Ear impression silicones must not be used for patients suffering from a perforated eardrum or injuries to the ear canal.

**Target Patient Group**

Patients who need ear care.

**Intended users**

ENT specialist, hearing aid audiologist.

**Product description**

The ear impression material consists of addition curing silicone with a soft, thixotropic consistency for impression technology which creates slight pressure. The material is practically shrink-proof and is physiologically well tolerated thanks to the platinum hardener system.

**Processing**

The ear canal and the eardrum must be examined before the impression is taken. No impression may be taken in the event of abnormal findings such as inflammation or perforation of the eardrum. Cleaning of the ear canal and removal of hairs by an ENT doctor is required in the event of significant earwax deposits and hairs in the ear canal. An impression pad must be placed at the end of the outer ear canal before the eardrum before the impression is taken in order to protect the eardrum. The two components of the impression material are removed from the containers in equal parts with measuring spoons and pressed onto a mixing tray. The two components are now rapidly kneaded together with a spatula until the mixture is a uniform color. For hygiene reasons, kneading by hand is only recommended with gloves (no latex gloves!). The kneaded mass is now filled into an impression syringe. The impression material is slowly applied into the ear canal and the concha. Please take care to ensure that the tip of the impression syringe is always in direct contact with the impression mass during application. Once the material has hardened, the impression

can be carefully removed from the ear. The ear canal must be examined again after the impression has been taken. Carefully close the tins after use.

### **Disinfection**

The impressions can be disinfected with commercially available disinfectants according to the manufacturer's instructions.

### **Caution**

- The product is exclusively intended for the above-mentioned range of use and shall only be used by professionally qualified and instructed persons.
- The setting times can be increased by lower temperatures or reduced by higher temperatures.
- Wear protective clothing! Set material is chemically resistant. Marks on your clothing cannot be removed.
- The chemical reaction between the two components may be inhibited by earwax and the use of latex gloves, hand lotions, cleaning agents, etc. for this reason, we recommend that you wear commercially available polyethylene gloves.
- Carefully close the tins after use.

### **Note for one-time use**

After taking an impression, the cured ear impression silicone cannot be used for taking another impression owing to the physical properties.

### **Disposal**

Disposal of contents/container in accordance with local/regional/national/international regulations and safety data sheet.

### **Hazard identification**

Hazard and safety notes can be found on the product label and the corresponding safety data sheet.

### **Serious incidents**

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

### **Batch number/use-by date**

The batch number and the use-by date are located on both the outer packaging as well as on each tub. Both components A and B of a packaging must only be used in this combination. When submitting a complaint about a product, please always quote the batch number. Do not use the product after expiration of the use-by date.

### **Technical data see page 2.**

**Innovation MediTech GmbH is not liable for any damages caused by improper application of the product.**

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egger Otoplastik + Labortechnik GmbH