egger flex / AB

60 shore

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

Instructions for use

Specifically developed for the 3D casting technique. Optimal air displacement with a low viscosity setting, no air bubble formation. Especially easy to demold, no adhesion to the cast mold. Permanently elastic, marked elastic recovery, high tear propagation resistance as well as elongation at break. High wear comfort though soft-flexible formulation, long-term stability, skin friendly & biocompatible. Medical device Class IIa.

1. Mixing and dosing

Extrude material from the cartridge with inject pistol S50. Place cartridge into inject pistol and remove closure cap by rotation. Eject a small amount of material (bleeding) until equal flow is being extruded from both orifices. Mount the mixing cannula by using the guides on the cannula on the cartridge and lock in place by rotating in the opposite direction. Feed material applying even pressure. Prior to use, dispense a small amount of material and visually check that the base and catalyst paste are mixed homogeneously. Only then dosing can be performed individually. Store cartridge with mixing cannula attached for subsequent use.

2. Application in the lab – manual

Prepare the ear impression for further processing and scan. Process the data for creation of the negative molds using appropriate software. Plaster or resins can be used as materials for the negative molds. Generally, plastic negative molds do not need any isolation, however, demolding can be simplified by rinsing with the separating agent Cast Separator e3D.CS01, followed by drying. Fill egger flex/AB into the negative mold without bubbles using the inject pistol (Fig. 1). Vulcanisation of the bubble-free injected material takes place at room temperature. To avoid air bubbles, it is recommended to harden under pressure (pressure vessel). Following vulcanisation, split or break up the negative form and remove the blank.

3. Final trimming and varnishing

For molding and surface trimming use the special egger EF-cutters (for silicones) and trimmer (Fig. 2). For smoothing the surface, use the air-drying egger silicone lacquers, Lacquer L (REF 37500 ff) and Lacquer L nano (with antibacterial effect, REF 37510 ff), Lacquer opaco (REF 37527) or the hot-vulcanizing Lacquer signo (REF 39412 ff) according to their instructions for use (Fig. 3).

Important working hints

- Do not combine with condensation curing silicones.
- Cured silicone materials are chemically inert – spots on clothing should be avoided.
- Latex gloves and latex contaminated surfaces as well as cerumen, cremes and resins may inhibit the setting reaction of egger flex/AB (we recommend standard gloves made of nitrile rubber or polyethylene).
- The cartridges are filled in a vacuum, during quality control thoroughly tested and only released in air bubble free conditions.
- Since under certain environmental conditions (e.g. temperature, air pressure) air bubbles may be produced afterwards in the cartridge, curing of the material in a pressure-container is generally recommended.
- Where movement is impeded or blocked, discard the cartridge. Do not apply force.

Safety advice

egger is not liable for any damage caused by improper application of the impression material.

Further information

Silicone based materials are proven a million times. On condition of a proper application, undesired effects are not to be expected. However, reactions of the immune system like allergies, irritations, cannot be absolutely excluded. In case of doubt, we recommend to make an allergy test before the application of the material.

For use by trained specialists.

Indications for use

Indirect fabrication of earmolds for connection to a BTE hearing system and for the fabrication of hearing protection earmolds.

Technical Data

Mixed volume: 50 ml (S50-cartridge)

Mixing ratio: 1:1

Product colors: colorless-transparent, reddish-transparent, green, pink-fluorescent, yellow-fluorescent

Mixing time: omitted (Automix system)

Working time: approx. 2 min. 30 sec. (from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity. Increased temperatures accelerate, decreased temperatures retard a.m. times.)

Setting time: approx. 15 min. (from beginning of mixing at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity. Increased temperatures accelerate, decreased temperatures retard a.m. times.)

Final hardness: 60 shore A

Recovery from deformation: > 99.8%

Linear dimensional change: 0.2%

Application: at 23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5% rel. humidity

Storage:

CAUTION:

Federal U.S. law restricts this device to sale by or on order of a licensed healthcare practitioner (or trained specialist personnel).

Please follow the product information and the material safety data sheet.