e3D.S01

beige, red plus, blue plus,
white, black

Allgemeine Produktinformation

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EN: Contains: aliphatic urethane acrylate, tripropyleneglycoldiacrylate, tetrahydrofurfuryl methacrylate THFMA purified grade, diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide, phenyl bis(2,4,6-trimethylbenzoyl)-phosphine oxide

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EN: Contains: Aliphatic urethane acrylate, tripropyleneglycoldiacrylate tetrahydrofurfuryl methacrylate THFMA purified grade, diphenyl(2,4,6-trimethylbenzoyl)phosphine oxide

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EN: May cause an allergic skin reaction. Causes serious eye irritation. Harmful to aquatic life with long lasting effects.

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EN: Avoid breathing dust/fume/gas/mist/vapours/spray. Wear protective gloves/protective clothing/eye protection/face protection. If on skin: Wash with plenty of water. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash before reuse. Dispose of contents/container according to local and national require

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Please follow the product information and
the material safety data sheet.

79612405-2023-08

General product information

e3D.S01 is a light-curing resin used for the generative, automatic fabrication of ITE shells based on systems with a light source of UV/LED 378–388 nm. e3D.S01 is based on (meth) acrylate and free of bisphenol A.

Packaging unit 1000 g

Bottle	beige	REF 43027
	red plus	REF 43026
	blue plus	REF 43025
	white	REF 43028
	black	REF 43029

Important note

e3D.S01 is declared as medical device and is to be used by trained specialist personnel only.

Processing

The light unit must be in an orderly condition and the generatively fabricated earmold/hearing protector must be completely cured. The polymerization product and its properties depend in particular, among other factors, on post-processing. Correct post curing is essential for biocompatibility. Homogenize/roll the material prior to processing, i.e. with a roller mixer.

Maximum curing depth at direct post-exposure:

- 1.1 mm for e3D.S01 black
- 1.3 mm for e3D.S01 beige
- 2 mm for e3D.S01 red plus and e3D.S01 white
- 3 mm for e3D.S01 blue plus

Please consider that the final color shade is achieved after a few days.

Storage

Polymerization can already be induced through a minimal influence of light. e3D.S01 is to be stored dry and protected from light. The container must always be kept tightly sealed and immediately be closed carefully after each use.

Due to too cold temperatures during storage, the material may go cloudy. Short heating up to 35 °C can reverse this occurrence without loss of quality.

Safety

- Please follow the instructions on the safety data sheet!
- Irritating to eyes and skin (sensitization is possible).
- Avoid direct contact with the liquid material and the components prior to post-curing.
- Be sure to use personal protective equipment (protective gloves and protective glasses) during processing.
- After contact with eyes rinse thoroughly with water immediately and consult a doctor.
- After contact with skin wash immediately with water and soap.

Adverse effects

Product may cause allergic reactions.

Contraindication

Contains (meth)acrylics and phosphinoxides! The ingredients may cause allergic reactions in predisposed persons. In such cases refrain from using the product.

Note

e3D.S01 is intended and optimized to be used in the hearing aid sector. It must be used in accordance with the instructions for processing and safety. Egger will not be responsible for damages caused by faulty or improper use of system and materials.

Instructions for processing

1. Modeling

Please consider the respective instructions of the CAD software manufacturer.

2. Construction process

The print job has to be executed according to the instructions for use of the respective system.

3. Post-processing

If possible commence post-processing immediately following the completed construction process. Recommendation: drip time of 10 minutes after the print job is completed.

4. Pre-cleaning

Prepare a suitable vessel with isopropyl alcohol (purity ≥ 98%). As soon as you have carefully removed the construction components from the platform, and have placed them in the vessel with isopropyl alcohol, put the vessel in an ultrasonic bath for 3 minutes. Afterwards, remove the rest of the support structure. You can achieve good results when using compressed air for the cleaning of the component openings, e.g. sound canals and additional bore holes.

5. Main cleaning process

Keep another vessel with fresh isopropyl alcohol (purity ≥ 98%) ready. Place the construction components in the vessel and put it into the ultrasonic bath for another 3 minutes. Afterwards, check the openings and additional bore holes for residues.

6. Drying

In the oven at 40 °C (104 °F) for 30 minutes to vaporize the solvent residues from the cleaning process.

7. Post curing

For post curing use the xenon photoflash unit Otofash G171 and execute post-processing under inert gas conditions with 2 x 2000 flashes (rotate components after 2000 flashes!). Expose the color »black« with 3 x 2000 flashes.

8. Surface finishing

Recommendation: Varnishing the surface with Egger LP/H lacquer. Optional: Polishing the surface mechanically or drum-polishing; in this case, the biocompatibility of the earmold/hearing protector might be affected.

Machine suitability

- evolupt SolFlex
- evolupt Microlay
- evolupt MiiCraft
- evolupt Asiga
- Systems with a light source of UV/LED 378–388 nm

CAUTION: Federal U.S. law restricts this device to sale by or on the order of a licensed healthcare practitioner (or trained specialist personnel).