egger EL 2

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REF 51210 (115 V)

egger EL 2

Please read the operating instructions carefully before working with the device! Please keep the operating instructions in a safe place for future reference!

Index	page
1 General	13
2 Safety Regulations	13
2.1 For a safety conscious work environment	13
3 Installation, Handling and Initial Operation	14
3.1 Scope of delivery	14
3.2 Technical data	14
3.3 Spare parts and accessories	14
3.4 Assembling and connecting the unit	14
3.5 Initial operation and handling	15
3.6 Programming the timer	15
4 Working with UV-polymerizates	15
4.1 Storage	16
4.2 Specific processing notes	16
4.3 Polymerization of LP materials	16
4.4 After-polymerization with LP cleaner	16
4.5 RLacquering with LP/H lacquer or LP/H lacquer antibac	16
4.6 Lacquering with LP/W lacquer	17
4.7 Modifications/repairs with modeling pastes	17
5 Maintenance	18
5.1 Exchanging light tubes	18
5.2 Work performance of the light tubes	18
5.3 Arrangement of the light tubes	18
5.3 Exchanging of reflecting foil	18
6 Warranty/Service/After-Sales Service	19

Key

Important Information

Caution

1 General

The light unit egger EL 2 was developed to polymerize light-curing resins for the manufacturing and surface lacquering of BTE and ITE earmolds. This unit can process all kinds of egger LP materials as well as materials of different consistency and pigmentation from another manufacturer.

The light tubes (UV and blue light tubes) are arranged on different levels and can be switched on and off depending on their use. Consequently, this guarantees an optimal surface and depths polymerization.

This device can simultaneously polymerize 6 ITE or 4 BTE earmolds when using negative molds measuring a diameter of 50 mm.

The light unit egger EL 2 was especially developed for the a. m. applications and is characterized by

- maintenance-free working
- fast exchanging of the UV tubes
- easy handling
- simple cleaning of the powder-coated metal casing

In order to complete the entire processing method, we can offer you a variety of devices and useful accessories out of our wide range of products. Request our catalog or visit our homepage www.egger-labor.com. In case you need personal advice, we will be at your disposal.

2 Safety Regulations

 The electrical light-curing unit egger EL 2 is used for the polymerization of UV-curing LP materials, modeling pastes and lacquers. A use other than described in the operating instructions may lead to injuries of the operator or of third persons. Furthermore, the device or other material assets may be damaged. Therefore, the device must be used in technically perfect condition only and the safety regulations must be strictly complied with.

2.1 For a safety conscious work environment

- Protect the unit from dust. Use the unit on dry working places only.
- Qualified and trained persons must use the device only.
- Before putting the device into operation, the user has to read the operating instructions and ask the manufacturer in case of any questions.
- Use the unit in perfect and clean condition only and with fully operational UVA blue fluorescent light tubes.
- Before switching on the device, the light protecting cover has to be attached.
- Do not cover the ventilator of the unit. Do not put any items through the venting slot.
- Connect the device to permissible sockets only.
- Please adhere to the acknowledged safety regulations.
- This device is used for the processing of UV-curing materials in a spectral range of 360 - 390 nm. Please also adhere to the processing manual for each material to be worked with.
- Do not pull the mains plug out off the socket at the mains cable.
- Should the mains cable be damaged, it must be replaced at once. Defective mains cables must not be used.
- The device must be repaired by the manufacturer only.

egger EL 2

- The unit must be switched off during service and maintenance operations (cleaning, replacing of UVA/blue fluorescent light tubes or reflecting foil). Pull the mains plug before!
- When exchanging use the original egger light tubes and reflectors only.
- Do not use any dissolvent when cleaning the device. Do not immerse the unit in liquids.
- Children must be kept away from the unit.

! On self-accomplished repairs or if other spare parts are used, all guarantee and liability claims will become extinct!

3 Installation, Handling and Initial Operation

3.1 Scope of delivery

Please check whether the delivery is complete and whether all parts are free from transport damage. Immediately inform us or the distributor about missing parts or transport damage.

The UV-curing unit egger EL 2 consists of:

- 1 light protection lid
- 1 mains cable
- 1 operating and processing instructions and can be put into operation right away.

3.2 Technical data of the light-curing unit egger EL 2

Power supply: 230 V +/- 10% (or 115 V)

Mains frequency: AC 50 Hz

Connection: mains plug Schuko (German) (or US)

Power input: AC 80 VA

3.3 Spare parts and accessories for the light-curing unit egger EL 2

UVA light tube	REF 50103
Blue fluorescent light tube	REF 50104
Rectangle socket	REF 50107
Light protection lid	REF 51202
Set of reflecting foils	REF 51204
Plastic plate transparent	REF 51203
Timer	REF 51109

TIP: We recommend the service set, REF 50141, which consists of 2 UVA and blue fluorescent light tubes each and spare reflecting foils.

3.4 Assembling and connecting the unit

The electrical light-curing unit egger EL 2 corresponds to the state-of-the-art standards and to the acknowledged safety regulations, nevertheless dangers can arise. Please observe the following, before putting the unit into operation:

- Qualified and trained persons must use the device only.
- Please adhere to the safety regulations accepted.
- Use the unit in perfect and clean condition only and by observing the user instructions.
- Position the device horizontally and place it in its suitable working position.
- Remove the protection foils and the light protection cover off the workbench.
- Protect the unit against abrasion.
- Check, whether the unit is dry.
- Check whether the device is in technically perfect condition, especially the mains cable!
- Plug in the mains cable.

3.5 Initial operating and handling for the light-curing unit egger EL 2

Mains switch

- Place: reverse side of the unit.
- Function: During longer breaks switch off the unit.

Switch for light levels (A1/A2)

- Place: front panel of the unit.
- Function: The switch settings are marked with symbols. Depending on the switch setting (symbol), the light level is switched on, when starting the timer. The associated control lamp turns on.

Control lamps for the light levels

- Position: Front plate of the unit (in the switches for the light levels).
- Function: Indication of the light levels.

Acoustic signal

- Position: In the timer
- Function: Indicates the end of the polymerization process

3.6 Programming the timer

Adjusting the operating cycle

Switch on the device (main switch).

- Briefly press the button "PROG" (1 x for hour, 2 x for minute und 3 x for second), this is confirmed by an acoustic signal.
- The corresponding position on the display flashes. By using the button (+/-) you can set the desired value. Within the second and minute mode, the value can be changed in steps of 10 by pressing the (+/-) button longer.
- Confirm with the button PROG, then the time set will be shown on the display.

Starting and stopping the timer

- Starting and stopping the timer is done by pressing the button "START/STOP".
- After expiration of the polymerization time set, two signals will sound and on the display the bell symbol will be shown.
- Then, the timer automatically goes to its previous time setting.

Reset the timer to 0:00:00

- Simultaneously press the 3 buttons briefly.

Reset the timer to its start time

- After possibly stopping the preset polymerization time by pressing the "START/STOP" button, the timer can be reset to its previous time set by pressing the "RESET" button.

System reset

- Simultaneously press the 3 timer buttons for about 4 seconds.

4 Working with the UV-Curing Materials

4.1 Specific storage regulations

The materials are to be stored in a cool and dry place (at a maximum room temperature of 23° C / 73.4° F). Make certain that the containers are closed at any time, when not in use and close them immediately after use. Otherwise, there might be the danger of an unintentional chemical reaction by action of light (curing of the material). The expiry date, when stored adequately, is indicated on the product.

TIP: We recommend keeping already opened packages of lacquers or instant adhesives in the refrigerator.

egger EL 2

4.2 Specific processing notes

Please avoid direct skin contact with the materials. Before using them, protect your hands by wearing gloves or using special skin protection products.

The material is to be used by authorized and trained staff only. Processing of the lacquers should be done beneath an absorbing device (e.g. egger S/AB Timer REF 61400). Always use the lacquering brushes for one type of lacquer only. Please adhere to the operation instructions.

Details on the disposal can be taken from our safety data sheets that can be downloaded as PDF file from our website: www.egger-labor.com.

4.3 Polymerization of LP materials

- After filling the negative mold with LP material, for which you can use duplicating gel (REF 40900ff) or transparent duplicating silicone for (REF 41000ff), please put the negative form onto the workbench of the unit. Please adhere to the particular processing instructions for manufacturing ITE molds.
- Close the unit with the light protection lid.
- Set the required polymerization time at the timer (see processing instructions for the particular LP materials).
- By pressing the “START/STOP” button, the polymerization will be started.
- As soon as the polymerization is finished, a signal will sound. The unit will switch off automatically.
- Remove the light protection lid and take the negative mold from the workbench for further processing.

4.4 After polymerization with LP cleaner

The next step is to post-cure. This is necessary to remove the inhibition layer (sticky surface).

- Place the rectangular container filled with LP cleaner (REF 31300ff) onto the workbench of the unit.
- Insert the earmold into the filled rectangular container.
- Please make sure that the earmolds do not get in contact with each other.
- Close the unit with the light protection lid.
- Set the required polymerization time on the timer (see processing instructions for each specific LP material).
- By pressing the “START/STOP” button, the polymerization will be started.
- As soon as the polymerization is finished, a signal will sound. The unit will switch off automatically.
- Remove the light protection lid. Now, the earmolds can be taken out of the rectangular container with tweezers, bent (REF 26300).
- Then, the earmolds are rinsed under water and blown dry with compressed air.

4.5 Lacquering with LP/H lacquer and/or LP/H lacquer antibac

By using the LP/H lacquer (REF 3062ff) or LP/H lacquer antibac, hard earmolds can be protected due to an antibacterial coating with high gloss effect

- First the part of the earmold that is to be processed always has to be roughened with Corundum paper, grit 180, (REF 40704).
- Then, fix the BTE earmold onto a piece of nylon tube (REF 97200) and onto the lacquering holder. With an ITE earmold, the special Erkogum (REF 30911) can be attached to the opened battery valve.
- Then, clean and degrease the surface carefully by using the provided hand-spray pump with a surface cleaner (REF 31600ff).

- Apply the lacquer evenly. Alternatively, you can also immerse the earmold (but no ITE device with built-in electronics or face-plates!) into a dipping jar filled with lacquer (REF 31400 or 31401). Please shake off any lacquer remnants.
- Then, let the lacquer cure in the egger light unit for about 3 – 5 minutes. The earmold can then immediately be inserted.

4.6 Lacquering with LP/W lacquer

For the surface finish the light-curing soft earmolds (made of LP/W) are always lacquered with LP/W lacquer (REF 30700ff). Furthermore, all presently known hard acrylate-based otoplastic materials can be coated with LP/W lacquer to improve the fitting of the earmold (“non-skid coating”).

- First the part of the otoplastic that is to be processed always has to be roughened with Corundum paper, grit 180, (REF 40704).
- Then, fix the BTE earmold onto a piece of nylon tube (REF 97200) and onto the lacquering holder. With an ITE earmold, the special Erkogum can be attached at the opened battery valve.
- Then, clean and degrease the surface carefully by using the provided hand-spray pump with a surface cleaner (REF 31600ff).
- Apply evenly the lacquer and cure it for 5 – 7 minutes in the light unit.

After the polymerization, the otoplastic is post-cured as described in section 4.4.

NOTE:

Otoplastics for ITEs with built-in electronics can also be lacquered with LP/W lacquer.

Make certain that the otoplastic is not put into the LP-Cleaner bath. To remove the inhibition layer (sticky surface) use a surface cleaner moistened cloth or polymerize the otoplastic with gas.

4.7 Modifications/Repairs with Modeling Pastes

The LP/H modeling pastes can be universally used for various modifications and repairs with all kinds of hard UV-curing polymerizates or acrylates. These are available in the colors transparent, beige, red and blue.

- First of all, the part of the otoplastic that is to be processed has always to be roughened with Corundum paper, grit 180, (REF 40700).
- Then, clean and degrease the surface carefully by using the provided hand-spray pump with a surface cleaner (31600ff).
- Directly apply the modeling paste out of the tube by means of a modeling instrument (REF 31215) and smooth the transitions to the ground.
- Then, depending on the coating thickness applied, let it cure for about 4 – 6 minutes in the light unit egger EL 2.
- After the polymerization of the surface, wipe off the otoplastic by using a surface cleaner moistened cloth.
- If required, rework the otoplastic with a cutter.

5 Maintenance

! Before each maintenance, the device must be unplugged from the mains supply!

5.1 Exchanging light tubes

On each exchanging of the tubes (2 x UVA light tubes, REF 50103 and 2 x blue fluorescent light tubes, REF 50104) always exchange them all. Please make sure you always have some spare tubes in stock.

- Pull the tubes evenly backwards out of the socket. Do not cant!
- Simply put the new tubes into the socket until they are snapped in.
- The unit is now ready for operation and can be connected to the mains supply.
- Used light tubes are special waste and the legal provisions for their disposal have to be adhered to.

5.2. Work performance of the light tubes

Its work performance decisively depends on the switching on frequency and the mains voltage. The work performance averages approx. 1000 hours.


As the in-use time of the light tube per polymerization is pretty short, we recommend:

- exchanging the tubes each $\frac{1}{2}$ year, if the unit is in a daily use for about 8 hours
- exchanging the tubes each year, if the unit is in use for not more than 4 hours per day

Abrasion of the light tubes can be noticed on elongated polymerization times and a dark discoloration along the plug and socket area. Used light tubes are hazardous waste. Please observe the current legal regulations.

5.3 Exchanging of reflecting foil

Please ensure the correct arrangement of the light tubes after exchanging these, as they negatively influence the polymerization results, when positioned incorrectly.

5.4 Exchanging of reflecting foil

The polymerization performance of the light unit depends on the cleanliness of the light chamber. The built-in reflectors (REF 51204) amplify the irradiation and always have to be clean, otherwise, the polymerization times can change. Therefore, exchange the foils.

6 Warranty/Service/After-Sales Service

The warranty period for our devices starts with the date of delivery and is 12 months. Excluded from the warranty are modifications of the unit such as lamp modules, reflectors as well as all other parts of the unit which are subject to wear and tear.

 In case of improper handling of the electric light-curing unit and/or turning motor all warranty and liability claims will become extinct!

We will be pleased to be of assistance to you in case of damage or other malfunctions:

egger

Otoplastik + Labortechnik GmbH

Aybühlweg 59

87439 Kempten/Germany

Phone international: +49 831 58113-60

Fax international: +49 831 58113-14

E-mail international: sales@egger-labor.de

Internet: www.egger-labor.com

Please indicate the following data:

Device type:

Electric light-curing unit egger EL 2
(REF 51200/230 V or REF 51210/115 V)

serial no.:

(see label)

egger EL 2