


egger flex/AB 25
egger flex/AB 40

Otoplastiksilikone
Earmould silicones

egger

DE	Gebrauchsanweisung	6	PT	Manual de instruções	32
EN	Instructions for use.....	8	RO	Instrucțiuni de utilizare.....	34
FR	Mode d'emploi.....	10	SV	Bruksanvisning	36
ES	Modo de empleo	12	SK	Návod na použitie.....	38
IT	Istruzioni per l'uso	14	SL	Navodila za uporabo.....	40
DA	Brugsanvisning	16	CS	Návod k použití	42
EL	Οδηγίες χρήσης.....	18	HU	Használati útmutató.....	44
HR	Upute za uporabu.....	20	RU	Инструкция по применению	46
LV	Lietošanas instrukcija.....	22	TR	Kullanım kılavuzu	48
LT	Naudojimo instrukcija	24	FI	Käyttöohje	50
NL	Gebruiksaanwijzing	26	ET	Kasutusjuhend	52
NO	Bruksanvisning	28	BG	Инструкция за употреба	54
PL	Instrukcja stosowania.....	30			

Technical Specifications

	Mixing time	Working Time	Total setting time	Final hardness
egger flex/AB 25	not applicable (Automix1 system)	approx. 3 min 30 sec.*	approx. 30 min.* at room temperature (approx. 23 °C / 73,4 °F)	approx. 25 Shore A
egger flex/AB 40	not applicable (Automix1 system)	approx. 3 min 30 sec.*	approx. 30 min.* at room temperature (approx. 23 °C / 73,4 °F)	approx. 40 Shore A

Processing:

At 23 °C ± 2 °C / 73,4 °F ± 35,6 °F,
50 ± 5% rel. humidity

Storage:

* from start of mixing at 23 °C ± 2 °C / 73,4 °F ± 35,6 °F,
50 ± 5% rel. humidity.
Higher temperatures shorten, lower temperatures extend
the indicated times.

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Earmould silicone

Indication

- Indirect fabrication of earmoulds for connection to a BTE hearing system
- Indirect fabrication of earplugs
- Serial production

Patient target group

Persons for whom an earmould is to be created.

Intended users

Hearing aid acoustician, earmould laboratory

Technical data

See page 5

1. Mixing and metering

The cartridges are emptied with the Automix1 gun (Fig. 1). After inserting the cartridge into the mixing gun, remove the cartridge closure by turning it. Extrude a small amount of material until silicone is being uniformly emitted from both outlet openings (Fig. 2). To attach the mixing cannula, pay attention to the guides on the cannula and cartridge. Lock by turning clockwise (Fig. 3). Emit material with uniform pressure. Spread a small amount of material before use and visually check that the base and catalyst pastes are homogeneously mixed (Fig. 4). Only then can the metering be carried out individually. After use, leave the mixing cannula as a closure on the cartridge until the next application.

2. Processing in the lab

3D-printing polymers, gel or plaster can be used as materials for the negative mould. 3D-printed negative mould generally do not require isolation, but demoulding can be simplified by rinsing with e3D.CS01 separating agent and subsequent drying. Fill egger flex/AB (indirect method) into the negative mould, taking care to avoid bubbles (Fig. 5). Vulcanization of the material injected without bubbles takes place at room temperature (Fig. 6). To avoid formation of air bubbles, curing under pressure (pressure pot) is recommended. After vulcanization, split or break the negative mould and withdraw the blank.

3. Final trimming and varnishing

Final trimming can be carried out with EF milling cutters, pumice stones, or abrasive sleeves. The egger silicone varnishes Lack L nano (with silver particles), Lack opaco (matting) or Lack signo (hot-vulcanizing) can be used as surface finishes in accordance with the respective processing instructions. (Fig. 7).

Safety information

- Only for the specified intended use by trained specialists.
- Do not bring into contact with condensation-curing silicones.
- Cross-linked silicon materials are chemically resistant – avoid stains on clothing.
- Latex gloves and latex-contaminated surfaces can inhibit curing (we recommend nitrile or polyethylene gloves).
- In case of stiffness or clogging, discard the cartridge; do not use force for application.
- The cartridges are filled under vacuum, thoroughly checked during quality inspection, and released only if free of air bubbles. Since under certain ambient conditions air bubbles may subsequently develop in the cartridge (e.g. temperature, air pressure), curing in the pressure pot is generally recommended.
- Wear suitable personal protective equipment during the finishing of the cured material.
- Refer to the relevant safety data sheet for hazard and safety information.

Notes

- DETAX does not accept liability for any damage caused by misuse of the impression material.
- Read and understand the safety data sheet!

For users and/or patients:

Any serious incidents occurring in connection with this product must be reported immediately to incident@detax.com and to the competent public authority of the Member State in which the user and/or patient resides.

Contraindication

Do not use this material in case of known allergies to any of the ingredients, or contact allergies.

Silicone materials have been tried and tested millions of times; undesired effects are not to be expected in case of proper application. However, as a matter of principle immune reactions such as allergies or irritations cannot be ruled out. In case of doubt, we recommend performing an allergy test before application.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

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egger Otoplastik + Labortechnik GmbH
DETAX GmbH