









**egger flex/AB 40 direct / indirect**  
**egger flex/AB 60 direct / indirect**

Otoplastiksilikon  
Earmould silicone

**egger**

<b>DE</b>	Gebrauchsanweisung .....	6	<b>PT</b>	Manual de instruções .....	32
<b>EN</b>	Instructions for use.....	8	<b>RO</b>	Instrucțiuni de utilizare.....	34
<b>FR</b>	Mode d'emploi.....	10	<b>SV</b>	Bruksanvisning .....	36
<b>ES</b>	Modo de empleo .....	12	<b>SK</b>	Návod na použitie.....	38
<b>IT</b>	Istruzioni per l'uso .....	14	<b>SL</b>	Navodila za uporabo.....	40
<b>DA</b>	Brugsanvisning .....	16	<b>CS</b>	Návod k použití .....	42
<b>EL</b>	Οδηγίες χρήσης.....	18	<b>HU</b>	Használati útmutató.....	44
<b>HR</b>	Upute za uporabu.....	20	<b>RU</b>	Инструкция по применению .....	46
<b>LV</b>	Lietošanas instrukcija.....	22	<b>TR</b>	Kullanım kılavuzu .....	48
<b>LT</b>	Naudojimo instrukcija .....	24	<b>FI</b>	Käyttöohje .....	50
<b>NL</b>	Gebruiksaanwijzing .....	26	<b>ET</b>	Kasutusjuhend .....	52
<b>NO</b>	Bruksanvisning .....	28	<b>BG</b>	Инструкция за употреба .....	54
<b>PL</b>	Instrukcja stosowania .....	30			


## Technical Specifications

	Mixing time	Working Time	Total setting time	Final hardness
<b>egger flex/AB 40 direct / indirect</b>	not applicable (Automix1 system)	approx. 1 min 15 sec*	<u>direct method:</u> approx. 7 minutes at body temperature	approx. 40 Shore A
			<u>indirect method:</u> approx. 15 – 20 minutes* at room temperature (approx. 23 °C)	
<b>egger flex/AB 60 direct / indirect</b>	not applicable (Automix1 system)	approx. 1 min 15 sec*	<u>direct method:</u> approx. 7 minutes at body temperature	approx. 60 Shore A
			<u>indirect method:</u> approx. 15 – 20 minutes* at room temperature (approx. 23 °C)	

### Processing:

At 23 °C ± 2 °C /  
73 °F ± 4 °F, 50 ± 5%  
rel. humidity

### Storage:

\* from beginning of mixing at  
23 °C ± 2 °C / 73 °F ± 4 °F, 50 ± 5 % rel. humidity.  
Higher temperatures shorten, lower temperatures extend the  
indicated times.

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Earmould silicone

## Indication

- Direct and indirect fabrication of earmoulds for connection to a BTE hearing system
- Manufacture of hearing protection earmoulds

## Patient Target Group

Persons for whom an earmould is to be created.

## Intended Users

ENT specialist, hearing aid acoustician, earmould laboratory

## Technical data

See page 5

### 1. Mixing and metering

The cartridges are emptied with the Automix1 pistol (Fig. 1). After inserting the cartridge into the mixing gun, remove the cartridge closure by turning it. Extrude a small amount of material until silicone is being uniformly emitted from both outlet openings (Fig. 2). To attach the mixing cannula, pay attention to the guides on the cannula and cartridge. Lock by turning clockwise (Fig. 3). Emit material with uniform pressure. Spread a small amount of material before use and visually check that the base and catalyst pastes are homogeneously mixed (Fig. 4). Only then can the metering be carried out individually. After use, leave the mixing cannula as a closure cap on the cartridge until the next application.

### 2. Processing with the direct method

Before taking the impression, carefully examine the auditory canal and the eardrum (Fig. 5). In the event of pathological findings (e.g. injuries, inflammations or perforation of the eardrum, tube), no impression may be taken. The ear canal must be cleaned and depilated. Insufficient cleaning may delay or inhibit the curing process due to the presence of cerumen. It is mandatory to protect the eardrum with an impression pad (e.g. blue secure) before taking an impression. Place the impression pad at the end of the outer ear canal in front of the eardrum (Fig. 6). egger flex/AB (direct/indirect method) can then be applied directly into the pretreated, tamponed ear, whereby care is to be taken to avoid bubbles (Fig. 7). During the injection process, make sure that the tip of the cannula remains at the level of the injected material. After curing, the impression can be carefully withdrawn (Fig. 8). Finally, the ear must be re-checked for intactness (Fig. 9). The final trimming and surface finishing are carried out as described in Section 4.

### 3. Indirect processing in the laboratory

3D-printing polymers, gel or plaster can be used as materials for the negative moulds. 3D-printed negative moulds generally do not require isolation. Fill egger flex/AB (direct/indirect method) into the negative mould, taking care to avoid bubbles (Fig. 10). Vulcanization of the material injected without bubbles takes place at room temperature. To avoid formation of air bubbles, curing under pressure (pressure pot) is recommended. After vulcanization, withdraw the blank from the negative mould.

#### 4. Final trimming and varnishing

Final trimming can be carried out with EF milling cutters, pumice stones, or abrasive sleeves. The egger silicone varnishes Lack L nano (with silver particles), Lack opaco (matting) or Lack signo (hot-vulcanizing) can be used as surface finishes in accordance with the respective processing instructions. (Fig. 11).

#### Safety information

- Only for the specified intended use by trained specialists.
- Do not bring into contact with condensation-curing silicones.
- Cross-linked silicon materials are chemically resistant – avoid stains on clothing.
- In case of stiffness or clogging, discard the cartridge; do not use force for application.
- Do not leave any residues of the silicone material in the ear.
- Latex gloves and latex-contaminated surfaces, cerumen, creams and polymeric materials can inhibit curing (we recommend nitrile or polyethylene gloves).
- The processing instructions and precautionary measures must be observed. Otherwise, irreparable damage to the hearing organ or eardrum may result.
- After each impression taking, careful examination of the auditory canal is necessary.
- The cartridges are filled under vacuum, thoroughly checked during quality inspection, and released only if free of air bubbles. Since under certain ambient conditions air bubbles may subsequently develop in the cartridge (e.g. temperature, air pressure), curing in the pressure pot is generally recommended.
- Wear suitable personal protective equipment during the finishing of the cured material.
- Refer to the relevant safety data sheet for hazard and safety information.

#### Notes

- DETAX does not accept liability for any damage caused by misuse of the impression material.
- Read and understand the safety data sheet!

#### For users and/or patients:

Any serious incident that has occurred in connection with this product must be reported immediately to [incident@detax.com](mailto:incident@detax.com) and to the competent public authority of the Member State in which the user and/or patient resides.

#### Contraindication

Do not use this material in case of known allergies to any of the ingredients, or contact allergies.

Silicone materials have been tried and tested millions of times; undesired effects are not to be expected in case of proper application. However, as a matter of principle immune reactions such as allergies or irritations cannot be ruled out. In case of doubt, we recommend performing an allergy test before application.

#### Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

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egger Otoplastik + Labortechnik GmbH
DETAX GmbH