LP/H Lacquer
LP/H Lacquer antibac
LP/W Lacquer

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Varnishing of hard earmoulds

Indication

Light-curing varnish for coating acrylic-based ITEs / hard earmoulds

Patient Target Group

Persons for whom an earmould is to be created.

Intended Users

Hearing aid acoustician, earmould laboratory

Curing times in the various light curing units

see page 2

General product information

LP/H Lack and LP/H Lack antibac

All currently known acrylic-based hard earmould materials can be coated with LP/H Lack and LP/H Lack antibac.

LP/W Lack

Currently, all hard acrylic-based earmould materials known to us can be coated with LP/W Lack to improve the fit of the earmould and seal the surface.

Polymerization of these light-curing varnishes can be carried out in the light units egger EL1, EL 1 plus N₂, EL 2, EL 3, EL 4 and eLED. LP3 with reversible motor, as well as with any equivalent light units.

Processing instructions

The varnishes may be applied only by trained staff and in accordance with the relevant safety data sheet. The varnishes must be handled under a fume hood. Disposal instructions can be found in the safety data sheet for this product.

Processing instructions

1. Preparation

Roughen the surface of the earmould, e.g. with corundum paper (grit 180) egger order no 40704.

PLEASE NOTE: To avoid inhomogeneities, you should always roughen and repaint the entire earmould surface.

Fix the earmould on a piece of egger nylon tubing, order no 97100 (for BTE earmoulds) or on a piece of Erkogum egger, order no 30911 (for ITE earmoulds).

To obtain a surface free from dust and grease, clean it with egger surface cleaner, order no 31600.

2. Handling

Apply the varnish with a varnishing brush, or dip the earmould and then shake off the excess.

3. Curing

Only one polymerization is necessary for the LP/H varnishes (see page 2).
With LP/W Lack, the surface must be re-treated after the first polymerization.

- a) After the first exposure, post-expose the earmoulds for 6 minutes in the LP-Cleaner bath.
- b) Rinse with water and dry with a paper towel.
- c) Then immerse the earmoulds for 3 minutes in a bath of isopropanol/water 1:1 (w/w) in an ultrasonic bath, dry with isopropanol, and dry at 37 °C for approx. 60 minutes.

Safety Information

- Only for the specified intended use by trained specialists.
- Avoid direct contact with the liquid material. Hands must be protected by wearing suitable gloves when processing the earmould and when handling the varnish.
- Wear safety goggles when working with the product!
- The varnishes must be handled under a fume hood.
- Refer to the relevant safety data sheet for hazard and safety information.

Notes

- DETAX shall not be held liable for any damage caused by misuse.
- Re-seal the container tightly after use.
- Protect from direct exposure to light during use.
- In order not to impair the quality of LP/H Lack, LP/H Lack antibac or LP/W Lack, do not fill varnish residues back into the original bottle after use.
- Read and understand the safety data sheet!

For users and/or patients:

Any serious incident that has occurred in connection with this product must be reported immediately to incident@detax.com and to the competent public authority of the Member State in which the user and/or patient resides.

Contraindication

Contains (meth)acrylates and phosphine oxide.

Ingredients of LP/H Lack, LP/H Lack antibac and LP/W Lack may cause allergic reactions in accordingly predisposed persons. In such a case, do not use the product any further.

Side Effects

Product may cause allergic reactions.

Disposal

Disposal of the contents/container must be carried out in accordance with the local/regional/national and international regulations.

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egger Otoplastik+ Labortechnik GmbH
DETAX GmbH