egger A/ITE vario
egger A/soft vario
egger A/rapidSWITCH vario
egger A/rapid PRO vario
egger A/II vario

Additionsvernetzendes
Ohrabformsilikon

Addition curing ear
impression silicone

egger

DE Gebrauchsanweisung	6	PT Instruções de uso	32
EN Directions for use	8	RO Instrucţiuni de utilizare	34
FR Mode d'emploi	10	SV Bruksanvisning	36
ES Instrucciones de uso	12	SK Návod na použitie	38
IT Istruzioni per l'uso	14	SL Navodila za uporabo	40
DA Brugsanvisning	16	CS Návod k použití	42
EL Οδηγίες χρήσης	18	HU Használati útmutató	44
HR Upute za uporabu	20	RU Указания к применению	46
LV Lietošanas instrukcija	22	TR Kullanım kılavuzu	48
LT Naudojimo instrukcija	24	FI Käyttöohje	50
NL Gebruiksaanwijzing	26	ET Kasutusjuhend	52
NO Bruksanvisning	28	BG Инструкция за употреба	54
PL Instrukcja użytkowania	30		

Technical Specifications

	Product color [base / catalyst]	Mixing volume	Consistency low high	Processing time
egger A/soft vario	blue/white	50 ml cartridge		≥ 30 sec
egger A/ITE vario	pink/white	50 ml cartridge		≥ 30 sec
egger A/II vario	red/light red	50 ml cartridge		≥ 30 sec
egger A/rapidPRO vario	orange/white	50 ml cartridge		≥ 30 sec
egger A/rapidSWITCH vario	lilac/white	50 ml cartridge		≥ 30 sec

The processing time refers to a room temperature of $23^{\circ}\text{C} \pm 1^{\circ}\text{C}$ ($73.4^{\circ}\text{F} \pm 1.8^{\circ}\text{F}$) and a relative air humidity of 50 %, the setting time to a body temperature of $37^{\circ}\text{C} \pm 1^{\circ}\text{C}$ ($98.6^{\circ}\text{F} \pm 1.8^{\circ}\text{F}$). Technical data according to DIN EN ISO 4823.

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Caution!

The impression of the outer ear canal may only be taken with impression material by persons with appropriate qualifications. The processing information and precautions which are set out in these directions for use must be observed under all circumstances. Failure to observe these may result in irreparable damage to the ear or the eardrum. If damage or undesirable side effects occur despite the warnings, we request that you report this with a detailed description of your procedure and the damage to or side effects for the patient.

Intended use

Addition curing ear impression silicone

Indications for use

Ear impression material for reproduction of the outer ear canal of a patient with hearing loss, Tinnitus aurium or impaired hearing for the provision and adjustment of a hearing device.

Contraindication

The product may not be used if the patient is allergic to one of the ingredients. There are no side-effects to be expected from this medical device if processed and used properly. However, immune reactions (e.g. allergies) or local disorders of sensation (e.g. irritation in the auditory canal) cannot be ruled out in principle. If you become aware of adverse side-effects – and in case of doubt – please be sure to advise us in any case and describe the circumstances and symptoms as accurately as possible. We will investigate any indication.

Patient groups

Ear impression silicones must not be used for patients suffering from a perforated eardrum or injuries to the ear canal.

Target Patient Group

Patients who need ear care.

Intended users

ENT specialist, hearing aid audiologist.

Product description

Addition curing ear impression silicone in a soft, stable/thixotropic consistency for impression technology which creates slight pressure (egger A/II vario) and for pressureless impression technology (egger A/soft vario, egger A/ITE vario, egger A/rapidPRO vario, egger A/rapidSWITCH vario). egger A/rapidSWITCH vario with color indicator indicates that the material has cured and the ear impression can be removed through a color change from purple to pink. egger A/rapidPRO vario and egger A/rapidSWITCH vario have particularly short setting times. The material is practically shrink-proof and is physiologically well tolerated thanks to the platinum hardener system.

Processing

The ear canal and the eardrum must be examined before the impression is taken. No impression may be taken in the event of abnormal findings such as inflammation or perforation of the eardrum. Cleaning of the ear canal and removal of hairs by an ENT doctor is required in the event of significant earwax deposits and hairs in the ear canal. An impression pad must be placed at the end of the outer ear canal before the eardrum before the impression is taken in order to protect the eardrum. The impression material from the double cartridge is automatically mixed at a 1:1 ratio. For this, place the double cartridge into the injector. Remove the cap and dispense a little material from the double cartridge in order to ensure that the two components come out evenly. Now place a mixing canula on the double cartridge. Check whether the material is dispensed evenly and mixed from the canula. To do this,

dispense a little material again before beginning taking the actual impression. The impression material is slowly applied into the ear canal and the concha. Please take care to ensure that the tip of the mixing canula is always in direct contact with the impression mass which has been dispensed. After applying the impression material, please leave the mixing canula on the double cartridge as a cap and replace it with a new mixing canula for the next impression. Once the material has cured, the impression can be carefully removed from the ear. The ear canal must be examined again after the impression has been taken.

Disinfection

The impressions can be disinfected with commercially available disinfectants according to the manufacturer's instructions.

Caution

- The product is exclusively intended for the above-mentioned range of use and shall only be used by professionally qualified and instructed persons.
- The setting times can be increased by lower temperatures or reduced by higher temperatures.
- Wear protective clothing! Set material is chemically resistant. Marks on your clothing cannot be removed.
- Dispense a little of the material before attaching the mixing canula.
- Never dispense into a mixing canula which has hardened.
- After use, leave the mixing canula on the double cartridge as a cap.
- The chemical reaction between the two components may be inhibited by earwax and the use of latex gloves, hand lotions, cleaning agents, etc. For this reason, we recommend that you wear commercially available polyethylene gloves.

Note for one-time use

After taking an impression, the cured ear impression silicone cannot be used for taking another impression owing to the physical properties.

Disposal

Disposal of contents/container in accordance with local/regional/national/international regulations and safety data sheet.

Hazard identification

Hazard and safety notes can be found on the product label and the corresponding safety data sheet.

Serious incidents

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

Batch number /use-by date

The batch number and the use-by date are located on both the outer packaging and the product. When submitting a complaint about a product, please always quote the batch number. Do not use the product after expiration of the use-by date.

Technical data see page 4/5.

Innovation MediTech GmbH is not liable for any damages caused by improper application of the product.

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egger Otoplastik + Labortechnik GmbH